**Summit Oxygen 'Everest' Silicone Extension Tube**

**Product Information:**

- Product Reference Code – ET12 ET 7
- Description – 3.5mm ID (6mm OD) extension tube
- Tubing length – 12 metres
- Material - 70 shore translucent silicone
- Non-sterile
- Autoclavable
- For use with medical grade oxygen only
- Kink-resistant internal profile

**Latex Free**

 **Caution:**

- Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury or death.

 **DO NOT SMOKE WHILST USING OXYGEN!**

- DO NOT USE near open flames or very hot ignition sources
- Do not store in direct sunlight.
- Do not store in high temperature locations.
- Store in dry conditions.
- Do not use if this packaging has been opened or damaged.

**Directions:**

- Attach one end of the oxygen tubing to the oxygen flow meter and the other to the oxygen cannula.
- Set the oxygen flow as prescribed and place the nasal prongs **briefly underwater** to see that oxygen is flowing properly.
- **If there are no bubbles appearing in the water, do not use the cannula, and recheck the tubing connections.**
- Once it is confirmed that oxygen is flowing, insert nasal tips into the nostrils.

**Hygiene: Product suitable for autoclaving. In absence of this facility follow the instructions below**

- To disinfect, place for 20 minutes in boiling distilled water at least every 2 weeks or more often.
- Dry by flowing oxygen through the tubing.
- Do not use the tubing if water is still in the tubing.
- To remove surface dirt, wipe the cannula with a damp cloth containing a very weak solution of detergent.

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