
DE Gebrauchsanweisung EN Instructions for Use FR Mode d'emploi
NL Gebruiksaanwijzing IT Istruzioni d'uso TR Kullanma Kilavuzu
ES Manual de instrucciones

LMT 65880a 05/2021 DE, EN, FR, IT, TR, ES, MX

CE 0197

**Löwenstein Medical
Technology GmbH + Co. KG**
Kronsaalweg 40
22525 Hamburg, Germany
T: +49 40 54702-0
F: +49 40 54702-461
www.loewensteinmedical.com

LENA, LENA NV









Full Face Mask

LMT 65880a

LÖWENSTEIN
medical

LÖWENSTEIN
medical

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¹ The mask materials deteriorate if exposed to aggressive detergents, for example. In individual cases, it may be necessary to replace mask parts sooner.

8.1 Materials


Only use the mask after consulting your physician if you are allergic to any of these substances:

Mask part	Material
Mask cushion, forehead cushion, O ₂ inlet, emergency exhalation valve	SI (silicone)
Mask body, forehead support, elbow (vented), valve retainer	PA (polyamide)
Retaining ring vented, retaining ring non-vented, rotating sleeve, headgear clip, ripcord clip (optional), adjusting element	POM (polyoxymethylene)
Elbow non-vented	PC (polycarbonate)
Headgear, ripcord (optional)	PA (polyamide), PU (polyurethane), PES (polyester)

No parts of the mask contain latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate).

9 Characteristic pressure/flow curve

The characteristic pressure/flow curve shows leakage flow as a function of therapy pressure.

10 Symbols on the packaging

Symbol	Description
	Lot number
	Order number
	Manufacturer and if necessary date of manufacture

	CE marking
	Permitted temperature range for transport and storage
	Keep out of sunlight
	Use by date
	Indicates the product is a medical device
	Unique device identifier
	Follow the instructions for use

11 Warranty

Löwenstein Medical Technology gives the customer a limited manufacturer warranty on new original Löwenstein Medical Technology products and any replacement part fitted by Löwenstein Medical Technology in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions are available on the website of the manufacturer. We will also send you the warranty conditions on request.

Please bear in mind that any claim to warranty and liability shall be void if neither the accessories recommended in the instructions for use nor genuine replacement parts are used.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Masks including accessories	6 months

12 Declaration of conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaaßweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of the Medical Device Regulations (EU) 2017/745. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

In the EU: As a user and/or patient, you must report any severe adverse events occurring in connection with the product to the manufacturer and to the responsible authorities.

1 Operation

The following figures show you how to put on, adjust, remove, dismantle and assemble the mask and also how to supply oxygen:

- 1 Putting on the mask
- 2 Adjusting the mask
- 3 Removing the mask
- 4 Dismantling the mask
- 5 Assembling the mask
- 6 Supplying oxygen

2 Introduction

2.1 Intended use

LENA is intended to be used as accessory by patients weighing more than 30 kg who have been prescribed non-invasive positive airway pressure (PAP) therapy such as CPAP/APAP, BiLevel or NIV. The mask is for single-patient re-use at home and for single and multi-patient re-use in the hospital/institutional environment.

2.2 Contraindications

The mask may not be used in the following situations: necessity for immediate intubation; loss of consciousness, acute vomiting.

The mask may be used in the following situations only with particular caution: pressure points and acute injuries to the skin of the face; skin allergies involving the face; deformities of the face or nasopharynx; acute pain affecting the face; cough reflex restricted or absent; claustrophobia; acute nausea.

If you are not certain whether one of these situations applies to you, please consult your attending physician or medical advisor. Please also observe the contraindications in the instructions for use of your therapy device.

2.3 Side effects

The following side effects may occur with use of the mask: nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing. Should such side effects occur, please contact your attending physician or medical advisor.

3 Safety

Risk of injury from re-inhaling CO₂!

If the mask is used incorrectly, CO₂ may be re-inhaled.

- ⇒ Do not close off the exhalation system of the mask.
- ⇒ Only put on the mask for an extended period if the therapy device is running.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Patients who are unable to take the mask off themselves must be monitored by qualified nursing staff.

Risk of injury if the mask slips!

If the mask slips or falls off, the therapy is ineffective.

- ⇒ Monitor patients with restricted spontaneous breathing.
- ⇒ Activate low pressure/leakage alarms on the therapy device.

Risk of fire from the use of oxygen in combination with flammable substances!

Supplying oxygen without a special safety device can lead to fire and injure people.

- ⇒ Always use an oxygen safety valve.
- ⇒ Follow the instructions for use for the oxygen supply system.
- ⇒ Set up oxygen sources at a distance of over 1 m from the device.
- ⇒ Do not smoke.
- ⇒ Do not use naked flame.
- ⇒ Ensure adequate ventilation.
- ⇒ Keep mask free from oil and grease.

Risk of injury from anesthetic gases and atomization of drugs!

Anesthetic gas and atomized drugs may escape through the exhalation valve and put third parties at risk.

- ⇒ Never use the mask during anesthesia.
- ⇒ Do not use the mask to atomize liquid drugs.

Risk of injury from lack of cleaning!

The mask may be contaminated which may put the patient at risk.

- ⇒ Clean the mask before using for the first time (see section entitled "Cleaning and hygiene treatment").
- ⇒ Clean the mask regularly.
- ⇒ When selecting the detergent, consider potential allergies.

Risk of injury due to mask being damaged!

Impact can result in parts of the mask becoming loose or breaking off; these may put the patient at risk.

- ⇒ Perform a visual inspection before every use and after every cleaning operation.

4 Product description

4.1 Overview

A diagram of the individual parts can be found on the title page.

1	Headgear	7	Elbow
2	Forehead cushion	8	Emergency exhalation valve (vented only)
3	Forehead support	9	Rotating sleeve (vented only)
4	Adjusting element	10	Mask body
5	O ₂ inlet	11	Headgear clip
6	Retaining ring	12	Mask cushion

4.2 Compatible devices

With many unit combinations, actual pressure in the mask will not correspond to prescribed therapy pressure even if the therapy device is displaying the correct pressure. Have the unit combination adjusted by a physician or specialist dealer so that actual pressure in the mask corresponds to therapy pressure. The therapy pressure required may vary between different mask types, so prescription of a suitable therapy pressure should in each case involve adjusting/adapting therapy to suit the mask type which is also going to be used during therapy itself.

4.3 Exhalation system

Vented variants

Vented variants (transparent elbow and rotating sleeve) have an integrated exhalation system. The retaining ring and mask body are shaped so that there is a gap between these parts. The exhaled air can escape through this gap.

Non-vented variants

Non-vented variants (blue retaining ring and elbow) do not have an exhalation system. Only use non-vented full-face masks in combination with therapy devices which have an active exhalation valve and have alarms and safety systems for any device failure. If you are using the non-vented variant with an external exhalation valve, please see the associated instructions for use.

4.4 Emergency exhalation valve

WARNING

Risk of asphyxia if emergency exhalation valve not working properly!

Residues may cause the valve to stick and lead to CO₂ being re-inhaled.
⇒ Check before every use that the openings of the emergency exhalation valve are clear.

If the device fails, the emergency exhalation valve opens so that the patient can breathe ambient air.

4.5 Ripcord (optional)

The ripcord obtainable as an option allows the mask to be released quickly and easily in emergency situations (see "SOS" illustration).

5 Cleaning and hygiene treatment

WARNING

Risk of injury from inadequate cleaning!

Residues may block the mask, impair the integrated exhalation system and jeopardize the success of therapy.
⇒ For patients with a compromised immune system or particular background of illness, disinfect mask parts daily following consultation with the physician.

5.1 Clean mask

1. Dismantle mask (see Figure 4).
2. Clean mask in accordance with the table below:

Action	Daily	Weekly
Wash mask parts in hot water and a mild detergent.	X	
Clean mask parts thoroughly with a cloth or a soft brush when washing. Or: Put the mask parts in the dishwasher, max. temperature 70 °C.		X
Wash headgear by hand.		X

3. Rinse all parts with clean water.
4. Allow all parts to air-dry.
5. Perform a visual inspection.
6. If necessary: Replace damaged parts.
7. Re-assemble mask (see Figure 5).

Discolored mask parts do not impair the functionality of the mask.

5.2 Hygiene treatment (clinical sphere)

Risk of infection from inadequate hygiene treatment!

On change of patient, inadequate hygiene treatment may lead to a risk of infection for the patient.

⇒ In the event of a change of patient, subject the mask to careful hygiene treatment in line with the "Hygiene treatment information" brochure. The brochure can be found on the manufacturer's website. We will also send you this brochure on request.

6 Disposal

You can dispose of all parts in domestic waste.

7 Troubleshooting

Fault	Cause	Remedy
Pain due to pressure on the face.	Mask too tight.	Loosen headgear slightly.
Draft in the eye.	Mask too loose.	Tighten headgear slightly.
	Mask does not fit.	Contact your specialist dealer.
Therapy pressure is not reached.	Mask not correctly adjusted.	Re-adjust mask (see Figure 2).
	Mask cushion damaged.	Replace mask cushion.
	Breathing circuit leaking.	Check connectors, check tubes properly located.
	Emergency exhalation valve defective.	Replace the mask.

8 Technical specifications

	vented	non-vented (NV)
Product class to directive MDR (EU) 2017/745	IIa	
Dimensions (H x W x D) Size S	155 mm x 100 mm x 95 mm	155 mm x 100 mm x 105 mm
Size M	165 mm x 100 mm x 95 mm	165 mm x 100 mm x 105mm
Size L	175 mm x 100 mm x 100 mm	175 mm x 100 mm x 110 mm
Weight Size S	137 g	135 g
Size M	141 g	139 g
Size L	150 g	148 g
Dead space Size S	246 ml	252 ml
Size M	288 ml	270 ml
Size L	326 ml	321 ml
Therapy pressure range	4 hPa - 35 hPa	
Tube connection: tapered connection to EN ISO 5356-1	Ø 22 mm (male)	Ø 22 mm (female)
Temperature range: Operation	+5 °C to +40 °C	
Transport and storage	-20 °C to +70 °C	
Flow resistance at 50 l/min	0.32 hPa	0.04 hPa
at 100 l/min	0.67 hPa	0.14 hPa
Flow resistance, anti-asphyxia valve Inspiration at 50 l/min:	0.6 hPa	-
Exhalation at 50 l/min:	0.8 hPa	-
Switching pressure Anti-asphyxia valve		
• Open:	0.5 hPa	-
• Close:	2.2 hPa	-
Quoted two-figure noise emission value to ISO 4871:		
- Sound pressure level	12 dB(A)	-
- Sound power level	20 dB(A)	-
- Uncertainty factor	3 dB(A)	-
Service life	5 years	
Useful life	up to 12 months ¹	
Standards applied	EN ISO 17510: 2020	

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