
FITTING YOUR MASK


1


2

3

4

F&P Eson™ 2 NASAL MASK


Use and Care Guide

F&P Eson™ 2

**Thank you for choosing
your F&P Eson 2.**

This nasal mask has been designed to make it easier for you to fit CPAP therapy into your life. This is your mask Use and Care Guide. Please read the entire guide before you first use this mask. Following the steps in this guide, and referring to the images and diagrams, will assist you with the assembly, fitting and daily use of your F&P Eson 2 Nasal Mask.

LANGUAGES

PAGE	LANGUAGE
1	English
11	Français (French)
21	Español (Spanish)
31	Português (Portuguese)
41	Deutsch (German)
51	Nederlands (Dutch)
61	Italiano (Italian)
71	Svenska (Swedish)
81	Suomi (Finnish)
91	Norsk (Norwegian)
101	Dansk (Danish)
111	תִּירָבֶּע (Hebrew)
121	العربية (Arabic)
131	Türkçe (Turkish)

Manufacturer
Fisher & Paykel Healthcare Ltd
15 Maurice Paykel Place
East Tamaki, Auckland 1013

PO Box 14 348, Panmure
Auckland 1741
New Zealand

Tel: +64 9 574 0100
Fax: +64 9 574 0158
Email: info@fphcare.co.nz
Web: www.fphcare.com

Australia
Fisher & Paykel Healthcare Pty Limited
36-40 New Street, PO Box 167
Ringwood, Melbourne
Victoria 3134, Australia

Tel: +61 3 9879 5022
Fax: +61 3 9879 5232

Austria
Tel: 0800 29 31 23
Fax: 0800 29 31 22

Benelux
Tel: +31 40 216 3555
Fax: +31 40 216 3554

China
Tel: +86 20 3205 3486
Fax: +86 20 3205 2132

Finland
Tel: +358 (0)405 406618
Fax: +46 (0)8 36 6310

France
Tel: +33 1 6446 5201
Fax: +33 1 6446 5221

Germany
Tel: +49 7181 98599 0
Fax: +49 7181 98599 66

India
Tel: +91 80 4284 4000
Fax: +91 80 4123 6044

Irish Republic
Tel: 1800 409 011

Italy
Tel: +39 06 7839 2939
Fax: +39 06 7814 7709

Japan
Tel: +81 3 5117 7110
Fax: +81 3 5117 7115

Korea
Tel: +82 2 6205 6900
Fax: +82 2 6309 6901

Northern Ireland
Tel: 0800 132 189

Russia
Tel and Fax: +7 495 782 21 50

Spain
Tel: +34 902 013 346
Fax: +34 902 013 379

Sweden
Tel: +46 8 564 76 680
Fax: +46 8 36 63 10

Switzerland
Tel: 0800 83 47 63
Fax: 0800 83 47 54

Taiwan
Tel: +886 2 8751 1739
Fax: +886 2 8751 5625

Turkey
Fisher Paykel Sağlık Ürünleri
Ticaret Limited Şirketi,
Alinteri Bulvari 1161/1 Sokak
No. 12-14, P.O. Box 06371 Ostim,
Ankara, Turkey
Tel: +90 312 354 34 12
Fax: +90 312 354 31 01

UK **EC** **REP**
Fisher & Paykel Healthcare Ltd
Unit 16, Cordwallis Park
Clivemont Road, Maidenhead
Berkshire SL6 7BU, UK

Tel: +44 1628 626 136
Fax: +44 1628 626 146

USA/Canada
Tel: +1 800 446 3908
or
+1 949 453 4000
Fax: +1 949 453 4001

INTENDED USE

The F&P Eson 2 Nasal Mask is intended to be used by individuals who have been diagnosed by a physician as requiring CPAP or Bi-Level therapy. The F&P Eson 2 Nasal Mask is intended for single patient adult (≥ 66 lb (30 kg)) use in the home and for multiple patient adult use in the hospital or other clinical setting where proper disinfection of the device can occur between patient uses.

OPERATING INSTRUCTIONS

- The operating pressure range of your mask is 4 to 30 cmH₂O.
- The operating temperature range of your mask is 5 to 40 °C (41 to 104 °F).
- Before using your mask each time, you should:
 - i. Inspect it for damage or deterioration. If there is any visible deterioration (cracking, tears, etc.), do not use your mask, and seek replacement part(s) from your healthcare provider.
 - ii. Inspect the exhaust holes on the Elbow **E**. Ensure air is flowing through these holes and do not use if blocked.
 - iii. Make sure path of air in mask is free from objects that may prevent breathing.

NOTE: Failure to follow the operating instructions above may compromise the performance and safety of the mask.

i MASK PARTS

Refer to the mask diagram on the inside cover of this guide.





- | | |
|-------------------------|-----------------------|
| A Silicone Seal | E Elbow |
| B Mask Frame | F Swivel |
| C Headgear | G Frame Socket |
| D Headgear Clips | H Diffuser |

Note: The Elbow **E** and the Headgear Clips **D** are not intended to be removed from the mask by the home user. If removed unintentionally, refer to technical information page for assembly.

ii FITTING YOUR MASK

1. Remove all packaging before using your mask.
2. Hold open the bottom straps of the Headgear **C** with the Mask Frame **B** facing downwards.
3. Guide the headgear **C** over your head so that the silicone seal **A** is resting on your forehead.
4. Place your hands on the bottom straps on either side of the mask frame **B** and pull out and downwards. Ensure the silicone seal **A** rests over your nose and the blue top strap is sitting on your forehead.
5. Gently tighten the headgear straps of the mask headgear **C** starting with the blue forehead straps followed by the bottom straps.

6. The forehead straps should rest horizontally on the forehead, just above the ears. If required the crown strap can be adjusted to achieve this.

Useful fitting tips:

- If you are experiencing mask leaks, you may need to lift the Silicone Seal **A** off your face and reposition it again.
- Ensure the mask's headgear straps are not fastened too tightly.
- The forehead straps should rest horizontally on the forehead, just above the ears.

DISASSEMBLY FOR CLEANING

1. Unhook the Headgear Clips **D** from the Mask Frame **B**.
2. Slide the forehead straps of the Headgear **C** downwards from the top of the Mask Frame **B**.
3. Remove the Silicone Seal **A** from the Mask Frame **B**.
4. Remove the Swivel **F** from the Elbow **E**.
5. Remove the Diffuser **H** from the Elbow **E**.

CLEANING YOUR MASK AT HOME

Wash your mask (excluding the Headgear **C** and Diffuser **H**) after each use:

1. Hand-wash your mask in soap dissolved in lukewarm water. Do not soak for more than 10 minutes.
2. Rinse thoroughly with fresh water, ensuring that all soap residue has been removed.
3. Leave all parts to dry out of direct sunlight before reassembling.

Wash your Headgear **C** and Diffuser **H** after 7 days of use:

1. Hand-wash your Headgear **C** in soap dissolved in lukewarm water. Do not soak your Headgear for more than 10 minutes.
2. Gently hand-wash the Diffuser **H** in soap dissolved in lukewarm water; avoid scrubbing the material. Do not soak your Diffuser for more than 10 minutes.
3. Rinse thoroughly with fresh water, ensuring that all soap residue has been removed.
4. Leave all parts to dry out of direct sunlight before reassembling.

Cleaning Cautions:

- Do not clean your mask with products containing alcohol, anti-bacterial agents, antiseptic, bleach, chlorine or moisturizer.
- Do not clean your mask in a dishwasher.
- Do not store your mask in direct sunlight.


The above actions may deteriorate or damage the mask and shorten its life. If there is any visible deterioration do not use your mask, but seek replacement part(s).

V


MASK ASSEMBLY

1. Push the Silicone Seal **A** onto the Mask Frame **B** until you hear it click firmly into place.
2. Push the Swivel **F** onto the Elbow **E**.
3. Place the Diffuser **H** onto the Elbow **E**.

1


2 3


To assemble headgear:

4. Thread the Headgear Clips **D** back onto Headgear **C**.
5. Thread the forehead straps of the Headgear **C** through the slots in the Mask Frame **B**.
6. Hook the Headgear Clips **D** onto the Mask Frame **B**.

4 5

6

OXYGEN/PRESSURE PORT CONNECTOR


If pressure readings and/or additional oxygen are required, an Oxygen/Pressure Port Connector is available (REF 900HC452). Please contact your healthcare provider.

TECHNICAL SPECIFICATIONS

F&P ESON 2 EXHAUST FLOW

The F&P Eson 2 Nasal Mask has a unique Diffuser System that is designed to diffuse air which is expelled from the mask. This reduces draft. The Diffuser does not significantly change the amount of air flowing from the mask; therefore, the mask can be used with or without the Diffuser.

Pressure (cmH ₂ O)	4	6	8	10	12	14	16	18	20	22	24	26	28	30
Flow with Diffuser (L/min)	19	25	29	34	38	41	45	48	52	55	58	61	63	66
Flow without Diffuser (L/min)	27	34	39	44	48	52	56	59	62	66	69	71	74	77

RESISTANCE TO FLOW

Pressure drop through mask with Diffuser @ 50 L/min: 0.19 ± 0.1 cmH₂O
 Pressure drop through mask with Diffuser @ 100 L/min: 0.80 ± 0.1 cmH₂O
 Pressure drop through mask without Diffuser @ 50 L/min: 0.17 ± 0.1 cmH₂O
 Pressure drop through mask without Diffuser @ 100 L/min: 0.73 ± 0.1 cmH₂O

MASK DEAD SPACE

Small: 69 cc Medium: 86 cc Large: 98 cc

SOUND

The sound power level of mask with Diffuser is 21.3 dBA, with uncertainty 2.5 dBA
 The sound pressure Level of the mask with Diffuser is 13.3 dBA, with uncertainty 2.5 dBA.
 The sound power level of mask without Diffuser is 31.4 dBA, with uncertainty 2.5 dBA
 The sound pressure Level of the mask without Diffuser is 23.4 dBA, with uncertainty 2.5 dBA.

The F&P Eson 2 Nasal Mask and accessories are not made with natural rubber latex.

MULTI-PATIENT CLEANING INSTRUCTIONS

Your F&P Eson 2 Nasal Mask is provided patient-ready.

To reuse between patients please follow the guidelines on the following website:
www.fphcare.com

If you do not have Internet access please contact your Fisher & Paykel Healthcare Representative.

HEADGEAR (C AND D)

1. Hand-wash the Headgear (C and D) in soap dissolved in lukewarm water. Air-dry out of direct sunlight.

DIFFUSER (H)

2. The Diffuser (H) is for single patient use. Replace between patients.

ALL OTHER PARTS

3. Clean using a mild alkaline detergent according to manufacturer's instructions e.g. soak in 1% Alconox for 10 minutes, gently scrub components in solution with a soft non-metallic brush, rinse in 5 liters of potable water, repeat rinse. Dry completely out of direct sunlight.

4. HIGH-LEVEL DISINFECTION

Pasteurize at
 75 °C for 30 mins,
 80 °C for 10 mins or
 90 °C for 1 min.

The F&P Eson 2 Nasal Mask has been validated to withstand 20 cycles of the above parameters. Failure to properly clean this device may result in inadequate disinfection. It is the user's responsibility to qualify any deviations from the recommended method of reprocessing.

ADVANCED DISASSEMBLY/ASSEMBLY

Please refer to the diagrams at the back of this guide.

1. **To disassemble Elbow (E) from Mask Frame (B):**
 - Position Elbow (E) so it is facing upward.
 - Pull Elbow (E) downward, away from Mask Frame (B), until Elbow (E) pops out.
2. **To assemble Elbow (E) back into Mask Frame (B):**
 - Position Elbow (E) so it is facing upward.
 - Push Elbow (E) at an angle and rotate into the Mask Frame (B) until the Elbow (E) pops in.
3. **To disassemble Headgear Clips (D) from the Headgear (C):**
 - Lie the flexed tabs flat against the bottom straps of the Headgear (C) and guide through the Headgear Clip opening.
 - Unhook Headgear Clips (D) from Mask Frame (B).
4. **To assemble Headgear Clips (D) with the Headgear (C):**
 - Hook Headgear Clips (D) back onto Mask Frame (B).
 - Thread the bottom straps of Headgear (C) through the Headgear Clips (D) with flexed tabs on Headgear (C) straps facing away from Mask Frame (B).

WARNINGS

Your mask should only be used with approved CPAP or Bi-Level equipment recommended by your physician or respiratory therapist. Your mask should not be used unless the CPAP or Bi-Level machine is turned on and operating properly. The exhaust holes associated with the mask should never be blocked.

- **EXPLANATION OF WARNING:** CPAP and Bi-Level machines are intended to be used with special masks, which have exhaust holes to allow continuous flow of air out of the mask.
 - i. Device switched on (with mask connected): When the CPAP or Bi-Level machine is turned on and functioning properly, new air from the CPAP or Bi-Level machine flushes the exhaled air out through the mask's exhaust holes.
 - ii. Device switched off (with mask connected): When the CPAP or Bi-Level machine is turned off, not enough fresh air will be provided through the mask, and exhaled air may be re-breathed. Re-breathing of exhaled air for longer than several minutes can lead to suffocation in some circumstances.
- At low CPAP or Bi-Level pressures the flow through the exhaust holes may be inadequate to clear all exhaled gas from the tubing. Some re-breathing of exhaled air may occur.
- Discontinue use and consult your healthcare provider or physician if the following occurs:
 - i. Mask discomfort or irritation
 - ii. Allergic reaction
- Do not use if vomiting or feeling nauseous.
- Use with oxygen: With your CPAP or Bi-Level machine, the oxygen flow must be turned off when the machine is not operating. Oxygen accumulated in your CPAP or Bi-Level machine enclosure will create a risk of fire. You should not smoke or be in the vicinity of someone smoking while using oxygen with your mask.
- Only use F&P Eson 2 Diffuser spare parts.
- Ensure the diffuser is completely dry before each use.

CAUTIONS

Use your mask only for its intended use as directed in this guide.

WARRANTY STATEMENT

Fisher & Paykel Healthcare warrants that the mask (excluding foams and material diffusers), when used in accordance with its instructions for use, shall be free from defects in workmanship and materials and will perform in accordance with Fisher & Paykel Healthcare's official published product specifications for a period of 90 days from the date of purchase by the end user. This warranty is subject to the limitations and exceptions set out in detail in www.fphcare.com/sleep-apnea/masks/warranty-mask/

DISPOSAL

Your mask does not contain any hazardous material and may be disposed of with general waste.

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